

Tynor Orthotics Pvt. Ltd.
Product Box Matter
Knee Cap Comfeel
Code: D 23
Size: S, M, L, XL, XXL



Tynor Knee Cap Comfeel is scientifically designed to provide compression and support to the knee. This aesthetically designed product improves blood circulation around the knee and provides warmth that helps relieve pain and inflammation.

Product Intro for Box only: Gives compression and support to the knee

Construction: This product is made in a state of art machine with a combination of Polyisoprene, Polyamide, Cotton and Elastane.

Applications:

- Knee Osteoarthritis
- Muscular weakness
- Swelling and effusion of the knee
- Post-surgical rehabilitation
- Post-trauma management
- Psychological support to geriatric patients

Features:

- Scientific design gives optimal compression
- Anatomical shape conforms to the knee anatomy
- Designer looks for improved aesthetics
- Airy knit design ensures breathability
- Stretchable spandex ensures a snug fit
- Skin-friendly nylon gives supreme comfort
- Soft patellar region reduces pressure on patella
- Sleek design - wearable even under clothes

Level of Support: Level 1

Instructions:

- Use this product under the guidance of a qualified doctor.
- Discontinue use and seek guidance of a qualified doctor, in case of
 - i. Impaired sensation
 - ii. Impaired blood flow
 - iii. Allergy rash or itching
 - iv. If pain increases or persists
- Consult your doctor in case of concurrent use with another medical device.
- Misuse, or improper use of the device can lead to adverse effects and decreased effectiveness
- Follow the instructions mentioned on the box carefully
- Keep out of reach of children
- Dispose in a responsible manner as per hospital guidelines/countries regulatory norms.

Product Care:

Symbols

- Hand wash with mild detergent and water below 30°C without wringing.

Compiled by:

Reviewed & approved by:

Ref. No :	Issue No : A	Rev No.: 00	Reviewed by	Approved by
N/733/02	Issue date : 08.09.16	Rev Date : NIL	M.R.	M.D.

- Dry in shade on a flat surface
- Do not dry clean.
- Do not Iron.
- Do not bleach
- Any serious incident in relation to this device is reported to the Tynor (contacts provided at back and bottom of box/manual) and the competent authority of the state.

EC Rep:

MT Promedt Consulting GmbH
Ernst-Heckel-Straße 7
66386 St. Ingbert
Germany
tel.: +49 (0) 6894 581020
fax: +49 (0) 6894 581021
e-mail: info@mt-procons.com
www.mt-procons.com

Bullet Points for front:

- Optimal compression
- Relieved patellar pressure
- Heat retention

Bullet Points for Back:

- Pleasing aesthetic
- Anti slip silicone strip
- Soft patellar region
- Pull-on design

Contents: Polyisoprene, Polyamide, Cotton, Elastane

SIZE	INCHES	CM
S	14.8 to 17.2	37.6 to 44
M	17.2 to 19.6	44 to 50
L	19.6 to 22	50 to 56
XL	22 to 24.4	56 to 62
XXL	24.4 to 26.8	62 to 68

Measure circumference around mid-thigh approximately 6 inches / 15.2 cm above the knee.

How to Use:

- With the broader end facing upwards, wear the sleeve with a simple pull-on action.
- Ensure that the loose patellar region comes at patella.

Compiled by:

Reviewed & approved by:

Ref. No :	Issue No : A	Rev No.: 00	Reviewed by	Approved by
N/733/02	Issue date : 08.09.16	Rev Date : NIL	M.R.	M.D.