


SPAN Seating Products Matrix

SPAN Seating Products Matrix		Treatment w/ Positioning	Skin Protection				General Use			
		Equalizer®	Ez Dish®	Geo-Matt® Contour	Geo-Matt® PRT	Geo-Matt®	Geo-Matt® Wedge	Gel-T®	Cross Cut	Geo-Wave®
Level of Care	General Use			✓			✓	✓	✓	✓
	Skin Protection	✓	✓	✓	✓	✓				
	Positioning Feature	✓	✓	✓			✓			✓
Core	Geo-Matt®			✓	✓	✓	✓			✓
	Minimizes Shearing			✓	✓	✓	✓			
	Gel Components	✓						✓		
Cover	Non-Slip Bottom		✓	✓	✓	✓	✓	✓	✓	✓
	Two-way Stretch Fabric Top	✓								
	Wipe Clean		✓	✓	✓	✓	✓	✓	✓	✓
General	Medicare Code *Bariatric Code	E2607	E2603	E2603	E2603 E2604*	E2603	E2601	E2601 E2602*	E2601	
	Therapeutic Weight Limit (lbs.)	250	250	250	250	250	250	250	250	300
	Bariatric Weight Limit (lbs.)	500			500			500		
	Warranty	2 yrs.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	18 mos.	1 yr.	18 mos.
	Seat and Back Cushion									✓
	Height of Cushion (Inches)	3.5"	3"	3"	3" standard 4" bariatric	4"	4" front 2" back	2.5"	3"	3"
	Viscoelastic Top	✓	✓							

Equalizer®



EZ Dish®



Geo Matt® Contour



Geo Matt® PRT




Geo Matt®



Geo Matt® Wedge



Gel-T®



Cross Cut



Geo-Wave®

