

In-Situ Repositioning Sheet

M921539

- Used for repositioning and lateral transfers
 - Used to help with effective pressure management, to reduce and prevent pressure sores
 - Size L - constructed from 3D Spacer fabric
 - Size XL - constructed from 3D Spacer fabric with COOLMAX® and COOLMAX freshFX® EcoMade
 - Can be used as part of the bed make-up when used in combination with a draw sheet
- Elastic on corners so that the sheet will stay in place on the mattress
 - Generous size making it suitable for all body dimensions
 - Wash temp 95°C / 203°F
 - Tumble dry low heat (max 50°C)



Product Codes and Specifications:




SWL
1000lbs



SWL
1200lbs

	3D Spacer	Silver Spacer with COOLMAX® and COOLMAX freshFX® EcoMade	A	B
L	M921539-L		56"	78"
XL		M921539-XL	66"	78"



Savaria, trading as Silvalea Ltd. Copyright © 2026. All rights reserved.